Stable Micro Systems

Stress

g
Kg
N
mN

g/mm²
Kg/cm²
N/m²
Kg/m²
Pa
mPa

DR

KPa
Pa.s

°C

dB

g

Ohm
mO

formulas

Balance AED PT100

Define

See The Feeling

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

The diagram on the left shows the hardware components of the texture analyzer, categorized into three main parts: 力臂 (Lever Arm), 传动装置 (Transmission Device), and 底座 (Base). To the right is a photograph of the physical instrument, a white and blue Stabbe Micro System texture analyzer.


Lotun Science Best testing & Analysis Corp.

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
The screenshot shows the software interface with a project list on the left and a graph of Force (g) vs. Time (s) on the right. The graph shows a single data point at approximately 0.5 seconds with a force of about 100g. Below the graph is a table with columns for Parameter, No., Value, and Units.

Parameter	No.	Value	Units


The top screenshot is titled "力量校正" (Force Calibration) and shows a dialog box for "Calibrate Force" with fields for Date, Measured Point, and Type. The bottom screenshot is titled "高度校正" (Height Calibration) and shows a graph with a yellow arrow pointing down, indicating the adjustment of the probe height.

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
Command	Program	Flags	Comment
Mark Value	1: Clear Graph Results		
Drop Anchor	2: Redraw		
Select Anchor	3: Search Forwards		
Search Forwards	4: Go to Minimum Time		
Search Backwards	5: Drop Anchor		
Go to Minimum Time	6: Set Threshold Force 5 g		
Go to Maximum Time	7: Go to Peak Top Force		
Set Threshold	8: F1 = Mark Value(Force (Curren... R		
Go to Peak Top	9: Drop Anchor		
Go to Peak Bottom	10: W1 = Area (Traditional)(Active...		
Area (Traditional)	11: Go to Force 5 g		
Gradient	12: Drop Anchor		
Mean	13: W2 = Area (Traditional)(Active...		
	14: R0 = W2 / W1	F	
	15: R = R0 * 100	F	
	16: Store Variable R as "弹性" "%"		
	17: Go to Absolute Maximum Valu...		
	18: F2 = Mark Value(Force (Curren... R		
	19: F = F2 - F1	F	
	20: Store Variable F as "韧性" "g"		
	21: Mark Value(Strain (%)) As 脆性 R		

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